

ARVO X3

Measurement		Filter					
Fluorescence	Excitation ※1	355/40nm	485/14nm				
	Emission ※2	460/25nm	535/25nm				
Absorbance		260/10nm	280/10nm	405/10nm	450/8nm	492/8nm	595/10nm
Luminescence		400~700nm		(BRET1) 460/25nm - 535/25nm		(BRET2) 410/80nm - 515/30nm	

※1 The filters (405/10nm, 450/8nm, 492/8nm and 595/10nm) are also installed but may require additional adjustments.

※2 The filters (410/80nm and 515/30nm) are also installed but may require additional adjustments.